

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.52	
Mean	9.41	Mean	13.42
SD	0.12	SD	0.2
CV	1.27	CV	1.52
Day-Run	Value	Day-Run	Value
3 - 1	9.6	3 - 1	13.3
4 - 1	9.5	4 - 1	13
5 - 1	9.3	5 - 1	13.6
6 - 1	9.3	6 - 1	13.1
8 - 1	9.3	8 - 1	13
9 - 1	9.5	9 - 1	13.6
10 - 1	9.3	10 - 1	13.4
12 - 1	9.5	12 - 1	13.5
13 - 1	9.4	13 - 1	13.4
15 - 1	9.4	15 - 1	13.5
16 - 1	9.3	16 - 1	13.4
17 - 1	9.2	17 - 1	13.4
20 - 1	9.5	20 - 1	13.5
22 - 1	9.6	22 - 1	13.7
23 - 1	9.5	23 - 1	13.5
24 - 1	9.4	24 - 1	13.4
25 - 1	9.3	25 - 1	13.4
26 - 1	9.5	26 - 1	13.4
27 - 1	9.4	27 - 1	13.6
29 - 1	9.3	29 - 1	13.3
30 - 1	9.6	30 - 1	13.8

Westgard	تفسیر	Calcium	آزمایش	اتوآنالایزر BS-800	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	9.42	SD	0.14	Exp. Date	1403/1/17	Lot. No	035

