

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	10.2	
2 - 1	10.52	
4 - 1	9.96	
5 - 1	9.34	
7 - 1	9.64	
8 - 1	10.2	
9 - 1	10.8	
11 - 1	9.5	
12 - 1	9.75	
13 - 1	10.85	
14 - 1	10.5	
16 - 1	8.92	
18 - 1	9	
19 - 1	9.11	
20 - 1	10.82	
21 - 1	10.5	
22 - 1	9.73	
26 - 1	8.95	
27 - 1	10.2	
28 - 1	10	
29 - 1	10.3	
30 - 1	9.8	

Westgard	تفسیر	Progesterone	آزمایش	SIEMENS Centaur CP	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	9.97	SD	0.6	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	9.94	SD	0.62	CV	6.2			

