

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	28.86
SD	0.62
CV	2.16
Day-Run	Value
1 - 1	27.4
3 - 1	29.1
4 - 1	27.7
6 - 1	28.7
7 - 1	28.6
8 - 1	29
10 - 1	30
11 - 1	29.3
13 - 1	28.7
14 - 1	29
15 - 1	28.3
16 - 1	28.9
20 - 1	29.5
21 - 1	29.1
22 - 1	29
23 - 1	28.6
24 - 1	28.2
28 - 1	29.2
29 - 1	29.1
30 - 1	29.7

High - CV : 1.54	
Mean	137.11
SD	2.11
CV	1.54
Day-Run	Value
1 - 1	136.8
3 - 1	140
4 - 1	135.6
6 - 1	137.4
7 - 1	136.3
8 - 1	136.2
10 - 1	133
11 - 1	136.7
13 - 1	137
14 - 1	133.7
15 - 1	139
16 - 1	138
20 - 1	140
21 - 1	138.7
22 - 1	134.1
23 - 1	139.3
24 - 1	137.5
28 - 1	135.8
30 - 1	140

Westgard	تفسیر	SGPT (ALT)	آزمایش	اتوآنالایزر BS-800	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	28.59	SD	0.76	Exp. Date	1403/1/17	Lot. No	035

