

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

| Day-Run | Value | Description |
|---------|-------|-------------|
| 1 - 1 | 85.6 | 1:2S |
| 3 - 1 | 81.7 | |
| 4 - 1 | 78.7 | |
| 6 - 1 | 80 | |
| 7 - 1 | 77.1 | |
| 8 - 1 | 77.9 | |
| 10 - 1 | 79.6 | |
| 11 - 1 | 81.6 | |
| 13 - 1 | 78 | |
| 14 - 1 | 81 | |
| 15 - 1 | 81.2 | |
| 16 - 1 | 81.2 | |
| 20 - 1 | 79.4 | |
| 21 - 1 | 79.5 | |
| 22 - 1 | 75.2 | 1:2S |
| 23 - 1 | 77.7 | |
| 24 - 1 | 76.7 | |
| 28 - 1 | 82 | 1:2S 4:1S |
| 29 - 1 | 78.3 | |
| 30 - 1 | 79 | |

| | | | | | | | | | |
|------------|--------|---------|--------|--------------------|--------|-----------|-----------|---------|-----|
| Westgard | تفسیر | Amylase | آزمایش | اتوآنالایزر BS-800 | دستگاه | 8 | ماه | 1402 | سال |
| Control | biorex | Mean | 80.15 | SD | 2.08 | Exp. Date | 1403/1/17 | Lot. No | 035 |
| Next Month | Mean | 79.57 | SD | 2.33 | CV | 2.93 | | | |

