

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	91.1	
2 - 1	90.3	
3 - 1	90.6	
4 - 1	91.8	
6 - 1	91.2	
7 - 1	91.7	
8 - 1	90.5	
9 - 1	89.9	1:2S
10 - 1	91.5	
10 - 2	91.9	
13 - 1	90.8	
14 - 1	91.1	
15 - 1	90.5	
16 - 1	91.4	
17 - 1	91	
20 - 1	91.06	
24 - 1	91.6	
25 - 1	91.1	
27 - 1	91.3	
28 - 1	91.5	
29 - 1	90.9	
30 - 1	90.5	

Westgard	تفسیر	C3	آزمایش	BS-480 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	90.93	SD	0.49	Exp. Date	1402/11/7	Lot. No	641-01-030/2
Next Month	Mean	91.06	SD	0.52	CV	0.57			

