

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	17.7	
2 - 1	17.9	
3 - 1	18.3	1:2S
4 - 1	18.1	
7 - 1	18.1	
8 - 1	17.9	
9 - 1	17.84	
11 - 1	18.06	
13 - 1	18.3	1:2S
14 - 1	17.8	
15 - 1	18	
16 - 1	17.6	1:2S
17 - 1	17.9	
20 - 1	17.8	
24 - 1	18.2	
22 - 1	17.9	
23 - 1	18.1	
24 - 2	17.8	
25 - 1	18	
27 - 1	17.9	
28 - 1	17.7	
29 - 1	17.8	
30 - 1	18	

Westgard	تفسیر	C4	آزمایش	BS-480 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	17.92	SD	0.14	Exp. Date	1402/9/14	Lot. No	641-01-030/2
Next Month	Mean	17.94	SD	0.19	CV	1.03			

