

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	50.9	
2 - 1	49.3	1:2S
3 - 1	48.5	
4 - 1	49.6	
6 - 1	47.9	
7 - 1	50.3	
8 - 1	50.3	
9 - 1	46.9	1:2S
10 - 1	50.2	
11 - 1	52.8	
13 - 1	51.8	
14 - 1	52.8	
15 - 1	47.4	
16 - 1	50.1	
20 - 1	53.4	
21 - 1	53.5	
22 - 1	53.4	
23 - 1	48.9	
24 - 1	50.9	
25 - 1	55	1:2S
27 - 1	47	
28 - 1	47.8	
29 - 1	48.3	

Westgard	تفسیر	CA 125	آزمایش	2000 ایمولایت XPI-	دستگاه	8	ماه	1402	سال
Control	Seronorm	Mean	50.49	SD	1.79	Exp. Date	1402/4/25	Lot. No	2011909
Next Month	Mean	50.3	SD	2.34	CV	4.65			

