

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.88	
Mean	127.55	Mean	218.18
SD	1.61	SD	4.1
CV	1.26	CV	1.88
Day-Run	Value	Day-Run	Value
1 - 1	125	1 - 1	211
3 - 1	126	3 - 1	216
4 - 1	127	4 - 1	221
6 - 1	129	6 - 1	223
7 - 1	128	7 - 1	221
8 - 1	130	8 - 1	218.6
10 - 1	126	10 - 1	213
11 - 1	130	11 - 1	220
13 - 1	129	13 - 1	220
14 - 1	128	14 - 1	222
15 - 1	126	15 - 1	213
16 - 1	128	16 - 1	223
20 - 1	128	20 - 1	215
21 - 1	127	21 - 1	220
22 - 1	125	22 - 1	211
23 - 1	129	23 - 1	224
24 - 1	126	24 - 1	216
28 - 1	129	28 - 1	220
29 - 1	129	29 - 1	221
30 - 1	126	30 - 1	215

Westgard	تفسیر	CPK	آزمایش	اتوآنالایزر BS-800	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	127.68	SD	1.89	Exp. Date	1403/1/17	Lot. No	035

