

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.68	
Mean	9.4	Mean	13.24
SD	0.12	SD	0.22
CV	1.28	CV	1.68
Day-Run	Value	Day-Run	Value
1 - 1	9.6	1 - 1	13.1
3 - 1	9.4	3 - 1	13.2
4 - 1	9.3	4 - 1	13.4
6 - 1	9.4	6 - 1	13.3
7 - 1	9.3	7 - 1	13.2
7 - 2	9.3	7 - 2	13.2
8 - 1	9.3	8 - 1	13
10 - 1	9.4	10 - 1	13.3
11 - 1	9.4	11 - 1	13.1
13 - 1	9.4	13 - 1	13.5
14 - 1	9.4	14 - 1	13
15 - 1	9.3	15 - 1	13.1
16 - 1	9.6	16 - 1	13.5
20 - 1	9.4	20 - 1	13
21 - 1	9.3	21 - 1	13.1
22 - 1	9.6	22 - 1	13.2
23 - 1	9.6	23 - 1	13.7
24 - 1	9.5	24 - 1	13.7
28 - 1	9.2	28 - 1	13
29 - 1	9.3	29 - 1	13
30 - 1	9.5	30 - 1	13.4

Westgard	تفسیر	Calcium	آزمایش	BS-800 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	9.41	SD	0.12	Exp. Date	1403/1/17	Lot. No	035

