

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	30.1	
2 - 1	29.5	
3 - 1	29.7	
4 - 1	29.7	
6 - 1	29.8	
7 - 2	30.5	
8 - 1	30.4	
9 - 1	28.7	
8 - 2	31.4	1:2S
11 - 1	30.4	
13 - 1	30.9	
14 - 1	28.6	
15 - 1	28.5	
16 - 1	30.1	
17 - 1	29.5	
20 - 1	28.9	
22 - 1	30.2	
23 - 1	30.4	
24 - 1	30.1	
25 - 1	28.5	
27 - 1	30.2	
28 - 1	29.5	
29 - 1	31	1:2S
30 - 1	29.8	

Westgard	تفسیر	Ceruloplasmin	آزمایش	BS-480 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	29.63	SD	0.64	Exp. Date	1406/9/14	Lot. No	641-01-030/2
Next Month	Mean	29.85	SD	0.79	CV	2.65			

