

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.83	
Mean	2.61	Mean	3.69
SD	0.04	SD	0.07
CV	1.55	CV	1.83
Day-Run	Value	Day-Run	Value
1 - 1	2.61	1 - 1	3.73
3 - 1	2.6	3 - 1	3.74
4 - 1	2.6	4 - 1	3.76
6 - 1	2.66	6 - 1	3.65
7 - 1	2.63	7 - 1	3.68
8 - 1	2.6	8 - 1	3.64
10 - 1	2.56	10 - 1	3.7
11 - 1	2.56	11 - 1	3.66
13 - 1	2.61	13 - 1	3.63
14 - 1	2.58	14 - 1	3.63
15 - 1	2.57	15 - 1	3.7
16 - 1	2.61	16 - 1	3.7
20 - 1	2.68	20 - 1	3.74
21 - 1	2.65	21 - 1	3.74
22 - 1	2.55	22 - 1	3.6
23 - 1	2.66	23 - 1	3.87
24 - 1	2.66	24 - 1	3.75
28 - 1	2.62	28 - 1	3.7
29 - 1	2.65	29 - 1	3.6
30 - 1	2.55	30 - 1	3.6

Westgard	تفسیر	Creatinine	آزمایش	اتوآنالایزر BS-800	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	2.58	SD	0.1	Exp. Date	1403/1/17	Lot. No	035

