

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
8 - 1	24.8	
9 - 1	25.2	
10 - 1	25.1	
11 - 1	23.8	
13 - 1	23.97	
14 - 1	23.8	
15 - 1	24.2	
17 - 1	25.3	
18 - 1	25.1	
20 - 1	25.3	
21 - 1	25.2	
23 - 1	25.8	
24 - 1	25.2	
25 - 1	22.8	
27 - 1	23.6	
28 - 1	25.2	
29 - 1	25.3	
1 - 1	25.6	
2 - 1	24.2	

Westgard	تفسیر	Cu,urine	آزمایش	Atomic Absorption Spectrometry	دستگاه	8	ماه	1402	سال
Control	آزمایشگاه	Mean	24.93	SD	1.18	Exp. Date	1402/8/30	Lot. No	
Next Month	Mean	24.71	SD	0.81	CV	3.29			

