

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.42	
Mean	73	Mean	227.9
SD	0.86	SD	3.23
CV	1.18	CV	1.42
Day-Run	Value	Day-Run	Value
1 - 1	74	1 - 1	227
3 - 1	71	3 - 1	230
4 - 1	73	4 - 1	233
6 - 1	73	6 - 1	230
7 - 1	72	7 - 1	231
8 - 1	73	8 - 1	226
10 - 1	73	10 - 1	230
11 - 1	74	11 - 1	225
13 - 1	75	13 - 1	225
14 - 1	73	14 - 1	227
15 - 1	73	15 - 1	233
16 - 1	74	16 - 1	226
20 - 1	73	20 - 1	225
21 - 1	73	21 - 1	226
22 - 1	73	22 - 1	232
23 - 1	72	23 - 1	230
24 - 1	72	24 - 1	225
28 - 1	73	28 - 1	226
29 - 1	73	29 - 1	230
30 - 1	73	30 - 1	221

Westgard	تفسیر	FBS	آزمایش	BS-800 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	73	SD	1.11	Exp. Date	1403/1/17	Lot. No	035

