

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	22.8	1:2S
2 - 1	23.7	
3 - 1	22.9	1:2S
4 - 1	24.6	
6 - 1	23.9	
7 - 1	26.3	
8 - 1	25.8	
9 - 1	23.4	
10 - 1	24.5	
11 - 1	25.1	
13 - 1	25.4	
14 - 1	25.6	
15 - 1	24.7	
16 - 1	23.5	
17 - 1	26.3	

Westgard	تفسیر	FSH	آزمایش	2000 ایمولایت XPI-	دستگاه	8	ماه	1402	سال
Control	Seronorm	Mean	25.17	SD	1.01	Exp. Date	1402/6/5	Lot. No	2011909
Next Month	Mean	24.57	SD	1.17	CV	4.77			

