

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	12.2	
2 - 1	12.5	
3 - 1	12.2	
3 - 2	12.6	
4 - 1	12.5	
5 - 1	12.1	
6 - 1	12.5	
7 - 1	12.7	
8 - 1	12.5	
9 - 1	12.2	
10 - 1	11.9	
11 - 1	12.3	
12 - 1	12.4	
13 - 1	12.4	
14 - 1	12.2	
15 - 1	12.7	
16 - 1	12	
17 - 1	12.9	1:2S
17 - 2	12.4	
18 - 1	12.01	
19 - 1	12.4	
20 - 1	12.4	
21 - 1	12.02	
22 - 1	12.6	
23 - 1	12.5	
24 - 1	12.1	
25 - 1	11.9	
26 - 1	11.95	
27 - 1	12.5	
28 - 1	12.9	1:2S
28 - 2	12.4	
29 - 1	12.3	
30 - 1	12.2	
30 - 2	12.5	

Westgard	تفسیر	HB	آزمایش	سل کانتیر Siemens Advia	دستگاه	8	ماه	1402	سال
Control	R&D	Mean	12.3	SD	0.27	Exp. Date	1402/8/25	Lot. No	
Next Month	Mean	12.35	SD	0.26	CV	2.14			

