

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	67	
3 - 1	68.4	
4 - 1	68.5	
6 - 1	67.4	
7 - 1	66.7	
8 - 1	70	1:2S
10 - 1	67.7	
11 - 1	68.5	
13 - 1	68.7	
14 - 1	68	
15 - 1	66.8	
16 - 1	67.1	
20 - 1	69.6	
21 - 1	70	
22 - 1	69.1	
23 - 1	68.3	
24 - 1	68.8	
28 - 1	67.1	
29 - 1	70	
30 - 1	69.2	

Westgard	تفسیر	HDL	آزمایش	اتوآنالایزر BS-800	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	68.3	SD	1.62	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	68.35	SD	1.1	CV	1.61			

