

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	135.1	
2 - 1	129.3	1:2S
3 - 1	140.2	
4 - 1	141.2	
6 - 1	139.6	
7 - 1	136.9	
8 - 1	144.2	
9 - 1	130.7	
10 - 1	133.5	
11 - 1	143.3	
13 - 1	142.1	
14 - 1	136.4	
15 - 1	140.5	
16 - 1	135.2	
17 - 1	136.5	
20 - 1	137.1	
22 - 1	141	
23 - 1	140	
24 - 1	145	1:2S
25 - 1	139.5	
27 - 1	139.3	
28 - 1	141.2	
29 - 1	140.8	
30 - 1	138.5	

Westgard	تفسیر	Ig A	آزمایش	BS-480 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	137.29	SD	3.75	Exp. Date	1402/11/7	Lot. No	641-01-030/2
Next Month	Mean	138.63	SD	3.92	CV	2.83			

