

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	742	1:2S
2 - 1	771	
3 - 1	765	
4 - 1	760.3	
6 - 1	763.9	
7 - 1	768.2	
8 - 1	783	1:2S
9 - 1	761.3	
10 - 1	759	
13 - 1	771.8	
13 - 2	770.6	
14 - 1	757	
15 - 1	764.2	
16 - 1	755.6	
17 - 1	772.2	
20 - 1	773	
22 - 1	750	
23 - 1	753	
25 - 1	750	
25 - 2	780.5	
27 - 1	760.5	
28 - 1	770.6	
29 - 1	750.9	
30 - 1	760.3	

Westgard	تفسیر	Ig G	آزمایش	BS-480 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	763.09	SD	8.9	Exp. Date	1402/11/7	Lot. No	641-01-030/2
Next Month	Mean	763.08	SD	10.08	CV	1.32			

