

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	45.2	
2 - 1	47.3	
3 - 1	55.4	
4 - 1	54.2	
6 - 1	54.2	
7 - 1	50.7	
8 - 2	51.5	
9 - 1	45.4	
10 - 1	60.3	1:2S
11 - 1	50.3	
13 - 1	54.2	
14 - 1	49.1	
15 - 1	45.2	
16 - 1	47.1	
17 - 1	49.5	
21 - 1	48.1	
20 - 1	48.1	
22 - 1	50.3	
23 - 1	55.3	
24 - 1	47.6	
25 - 1	56.1	
27 - 1	50.6	
28 - 1	49.3	
29 - 1	53.2	
30 - 1	55.2	

Westgard	تفسیر	Ig M	آزمایش	BS-480	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	49.26	SD	3.81	Exp. Date	1402/11/7	Lot. No	641-01-030/2
Next Month	Mean	50.94	SD	3.92	CV	7.7			

