

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	230	
3 - 1	230.7	
4 - 1	230	
6 - 1	234	
7 - 1	233.9	
8 - 1	235	
10 - 1	230	
11 - 1	231.5	
13 - 1	238.6	
14 - 1	231.1	
15 - 1	228.2	
16 - 1	231.4	
20 - 1	226.9	
21 - 1	230.4	
22 - 1	236.5	1:3S
23 - 1	230.6	
24 - 1	225.9	1:2S
28 - 1	238	1:2S 2:2S
29 - 1	233	
30 - 1	236	

Westgard	تفسیر	Iron	آزمایش	BS-800 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	234.01	SD	3.79	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	231.38	SD	3.65	CV	1.58			

