

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.54	
Mean	312.94	Mean	431.34
SD	8.52	SD	6.65
CV	2.72	CV	1.54
Day-Run	Value	Day-Run	Value
1 - 1	306	1 - 1	427
3 - 1	308	3 - 1	424
4 - 1	298	4 - 1	438
6 - 1	322	6 - 1	428
7 - 1	309	7 - 1	425
8 - 1	311	8 - 1	438.1
10 - 1	310	10 - 1	430
11 - 1	314.8	11 - 1	443.7
13 - 1	320	13 - 1	426
14 - 1	300	14 - 1	431
15 - 1	309	15 - 1	435
16 - 1	310	16 - 1	427
20 - 1	310	20 - 1	438
21 - 1	308	21 - 1	428
22 - 1	325	22 - 1	423
23 - 1	310	23 - 1	436
24 - 1	311	24 - 1	439
28 - 1	320	28 - 1	440
29 - 1	330	29 - 1	430
30 - 1	327	30 - 1	420

Westgard	تفسیر	LDH	آزمایش	اتوآنالایزر BS-800	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	309.54	SD	8.79	Exp. Date	1403/1/17	Lot. No	035

