

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	29.6	1:2S
2 - 1	32.3	
3 - 1	32.7	
4 - 1	32.1	
6 - 1	32	
7 - 1	34.7	
8 - 1	32.7	
9 - 1	32.1	
10 - 1	32	
11 - 1	34.2	
13 - 1	34.3	
14 - 1	30.5	1:2S
15 - 1	31.5	
16 - 1	31.4	
17 - 1	34.6	
20 - 1	31.8	
21 - 1	33.5	
22 - 1	33.5	
23 - 1	30.6	
24 - 1	32.3	
25 - 1	33.9	
27 - 1	34.4	
28 - 1	30	1:2S
29 - 1	30.9	
30 - 1	29.8	1:2S

Westgard	تفسیر	LH	آزمایش	2000 ایمولایت XPI-	دستگاه	8	ماه	1402	سال
Control	Seronorm	Mean	33	SD	1.23	Exp. Date	1403/5/31	Lot. No	2011909
Next Month	Mean	32.3	SD	1.55	CV	4.8			

