

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	43.6	
2 - 1	44.2	
3 - 1	42.9	
4 - 1	43.8	
6 - 1	42.4	1:2S
7 - 1	43.4	
8 - 1	43.6	
9 - 1	43.9	
10 - 1	42.9	
11 - 1	43.1	
13 - 1	43.2	
14 - 1	44.1	
15 - 1	43.5	
16 - 1	44.5	1:2S
17 - 1	43.5	
20 - 1	42.8	
22 - 1	42.3	1:3S
23 - 1	44.2	
24 - 1	43.2	
25 - 1	44.3	
27 - 1	43.5	
28 - 1	44	
29 - 1	44.8	1:3S
30 - 1	43.2	

Westgard	تفسیر	Lactate(Lactic Acid)	آزمایش	BS-480 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	43.43	SD	0.45	Exp. Date	1402/12/14	Lot. No	
Next Month	Mean	43.54	SD	0.55	CV	1.26			

