

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	32	
2 - 1	33	
3 - 1	34	
4 - 1	31	
6 - 1	33	
7 - 1	34	
8 - 1	32	
9 - 1	33	
10 - 1	35	1:2S
10 - 2	34	
11 - 1	32	
12 - 1	33	
13 - 1	31.5	
14 - 1	32	
15 - 1	33	
16 - 1	32	
17 - 1	33	
18 - 1	32	
19 - 1	33	
20 - 1	33	
21 - 1	34	
22 - 1	32	
23 - 1	31	
24 - 1	33	
25 - 1	33	
27 - 1	34	
28 - 1	32	
29 - 1	32	
30 - 1	32	

Westgard	تفسیر	MCH	آزمایش	سل کانتیر Siemens Advia	دستگاه	8	ماه	1402	سال
Control	R&D	Mean	32.2	SD	0.96	Exp. Date	1402/8/25	Lot. No	
Next Month	Mean	32.71	SD	0.98	CV	2.99			

