

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	86	
2 - 1	86.7	
3 - 1	86.5	
4 - 1	85.9	
5 - 1	85	
6 - 1	87	
7 - 1	87	
8 - 1	85	
9 - 1	86	
10 - 1	86.5	
11 - 1	86.8	
13 - 1	86.5	
14 - 1	87	
15 - 1	86	
16 - 1	87	
17 - 1	88	
18 - 1	86	
19 - 1	86	
20 - 1	86	
21 - 1	86.5	
22 - 1	87	
23 - 1	87.1	
24 - 1	85	
27 - 1	87	
28 - 1	86	
29 - 1	87	
30 - 1	87	

Westgard	تفسیر	MCV	آزمایش	سل کانتیر Siemens Advia	دستگاه	8	ماه	1402	سال
Control	R&D	Mean	86.56	SD	0.81	Exp. Date	1402/8/25	Lot. No	
Next Month	Mean	86.43	SD	0.72	CV	0.83			

