

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	31	
2 - 1	35	
3 - 1	36	
4 - 1	29	
6 - 1	33	
7 - 1	34	
8 - 1	36	
9 - 1	32	
10 - 1	31	
11 - 1	35	
13 - 1	33	
14 - 1	30	
15 - 1	33	
16 - 1	37	
17 - 1	36	
18 - 1	32	
21 - 1	32	
22 - 1	35	
23 - 1	33	
24 - 1	35	
25 - 1	30	
27 - 1	33	
28 - 1	35	
29 - 1	30	

Westgard	تفسیر	PLT	آزمایش	سل کانتیر Siemens Advia	دستگاه	8	ماه	1402	سال
Control	R&D	Mean	33	SD	2.98	Exp. Date	1402/8/25	Lot. No	
Next Month	Mean	33.17	SD	2.26	CV	6.81			

