

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 2.37	
Mean	3.43	Mean	7.98
SD	0.12	SD	0.19
CV	3.55	CV	2.37
Day-Run	Value	Day-Run	Value
1 - 1	3.6	1 - 1	7.8
3 - 1	3.5	3 - 1	8.2
4 - 1	3.7	4 - 1	8.1
6 - 1	3.3	6 - 1	8.2
7 - 1	3.4	7 - 1	7.8
8 - 1	3.6	8 - 1	7.9
10 - 1	3.3	10 - 1	7.8
11 - 1	3.5	11 - 1	7.8
13 - 1	3.3	13 - 1	8.3
14 - 1	3.4	14 - 1	7.9
15 - 1	3.4	15 - 1	8
16 - 1	3.4	16 - 1	8.1
20 - 1	3.6	20 - 1	7.8
21 - 1	3.4	21 - 1	8.2
22 - 1	3.3	22 - 1	7.6
23 - 1	3.3	23 - 1	8
24 - 1	3.4	24 - 1	7.8
28 - 1	3.5	28 - 1	8.2
29 - 1	3.3	29 - 1	8
30 - 1	3.4	30 - 1	8

Westgard	تفسیر	Phosphorus ( P )	آزمایش	اتوآنالایزر BS-800	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	3.47	SD	0.11	Exp. Date	1403/1/17	Lot. No	035

