

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	4.7	
2 - 1	4.8	
3 - 1	4.6	
4 - 1	4.7	
6 - 1	4.65	
7 - 1	4.69	
8 - 1	4.8	
9 - 1	4.67	
10 - 1	4.73	
11 - 1	4.81	
13 - 1	4.6	
14 - 1	4.7	
15 - 1	4.68	
16 - 1	4.65	
17 - 1	4.77	
20 - 1	4.7	
22 - 1	4.9	1:3S
23 - 1	4.72	1:3S
24 - 1	4.8	
25 - 1	4.7	
27 - 1	4.75	
28 - 1	4.7	

Westgard	تفسیر	Potassium (K)	آزمایش	AUDICOM	دستگاه	8	ماه	1402	سال
Control	زیست شیمی	Mean	4.7	SD	0.06	Exp. Date	1404/1/14	Lot. No	210202-17
Next Month	Mean	4.71	SD	0.06	CV	1.33			

