

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	17.2	
2 - 1	17.1	
3 - 1	18.8	1:3S
4 - 1	19.2	
6 - 1	19.5	
7 - 1	19.7	
8 - 1	18.4	
9 - 1	17.8	
10 - 1	17.5	
11 - 1	18.1	
13 - 1	18.9	
14 - 1	18.1	
15 - 1	17.9	
17 - 1	17.3	
20 - 1	19.4	
21 - 1	18	
22 - 1	17.1	
23 - 1	18.1	
25 - 1	17.8	
27 - 1	18.3	
28 - 1	17.7	
29 - 1	17.9	
30 - 1	17.8	

Westgard	تفسیر	Prolactin	آزمایش	2000 ایمولایت XPI-	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	18.42	SD	0.67	Exp. Date	1403/2/10	Lot. No	365-02-028
Next Month	Mean	18.16	SD	0.77	CV	4.22			

