

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	6.98
SD	0.08
CV	1.19
Day-Run	Value
1 - 1	6.9
3 - 1	6.9
4 - 1	6.9
6 - 1	6.9
7 - 1	7
8 - 1	7
10 - 1	6.8
11 - 1	7
13 - 1	7.1
14 - 1	7.1
15 - 1	6.9
16 - 1	7.1
20 - 1	7.1
21 - 1	7
22 - 1	7
23 - 1	7
24 - 1	6.9
28 - 1	7
29 - 1	7
30 - 1	6.9

High - CV : 1.58	
Mean	8.27
SD	0.13
CV	1.58
Day-Run	Value
1 - 1	8.2
3 - 1	8.3
4 - 1	8.4
6 - 1	8.3
7 - 1	8.2
8 - 1	8.1
10 - 1	8.5
11 - 1	8.4
13 - 1	8.3
14 - 1	8.3
15 - 1	8.3
16 - 1	8.4
20 - 1	8.1
21 - 1	8.2
22 - 1	8.2
23 - 1	8.4
24 - 1	8.2
28 - 1	8.1
29 - 1	8.4
30 - 1	8

Westgard	تفسیر	Protein Total	آزمایش	اتوآنالایزر BS-800	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	6.96	SD	0.07	Exp. Date	1403/1/17	Lot. No	035

