

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
3 - 1	97.6	
4 - 1	97.2	
6 - 1	98.1	
7 - 1	97.6	
8 - 1	98.2	
9 - 1	98.8	
11 - 1	97.9	
13 - 1	98.6	
14 - 1	98.5	
15 - 1	99.2	
16 - 1	98.2	
17 - 1	98.7	
20 - 1	97.5	
23 - 1	98.5	
24 - 1	99	
25 - 1	99.2	
27 - 1	99.2	
28 - 1	99.1	
29 - 1	98.6	
1 - 1	99.6	
2 - 1	96.3	

Westgard	تفسیر	Selenium	آزمایش	Atomic Absorption Spectrometry	دستگاه	8	ماه	1402	سال
Control	آزمایشگاه	Mean	98	SD	1.4	Exp. Date	1402/8/30	Lot. No	
Next Month	Mean	98.36	SD	0.8	CV	0.82			

