

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	148.2	
2 - 1	147.5	
3 - 1	146.7	
4 - 1	147.5	
6 - 1	147.3	
7 - 1	146.9	
8 - 1	147.2	
9 - 1	148.2	
10 - 2	147.3	
11 - 1	146.9	
13 - 1	147.3	
14 - 1	148.1	
15 - 1	146.9	
16 - 1	147	
17 - 1	148.1	
20 - 1	148	1:3S
22 - 1	148	
23 - 1	148.3	4:1S
24 - 1	147.3	
25 - 1	148.2	
27 - 1	147.5	
28 - 1	148.1	
29 - 1	146.8	
30 - 1	147.2	

Westgard	تفسیر	Sodium (Na)	آزمایش	AUDICOM	دستگاه	8	ماه	1402	سال
Control	زیست شیمی	Mean	147.14	SD	0.6	Exp. Date	1404/1/14	Lot. No	2103061-4
Next Month	Mean	147.5	SD	0.53	CV	0.36			

