

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	1.02
SD	0.07
CV	7.13
Day-Run	Value
1 - 1	1.03
3 - 1	0.9
4 - 1	0.98
4 - 1	0.98
6 - 1	1.04
7 - 1	1.1
8 - 1	1.12
10 - 1	1
11 - 1	1
13 - 1	1.1
14 - 1	0.92
15 - 1	1.05
16 - 1	0.99
20 - 1	1.05
21 - 1	1
22 - 1	1.1
23 - 1	0.89
24 - 1	1.07
28 - 1	1.12
29 - 1	0.9
30 - 1	1

High - CV : 2.39	
Mean	4.27
SD	0.1
CV	2.39
Day-Run	Value
1 - 1	4.23
3 - 1	4.11
4 - 1	4.5
6 - 1	4.38
7 - 1	4.3
8 - 1	4.3
10 - 1	4.37
11 - 1	4.15
13 - 1	4.3
14 - 1	4.15
15 - 1	4.3
16 - 1	4.15
20 - 1	4.35
21 - 1	4.3
22 - 1	4.3
23 - 1	4.13
24 - 1	4.3
28 - 1	4.37
29 - 1	4.2
30 - 1	4.2

Westgard	تفسیر	Bili T	آزمایش	BS-800 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	1	SD	0.07	Exp. Date	1403/1/17	Lot. No	035

