

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.69	
Mean	174.2	Mean	232.35
SD	2.35	SD	3.92
CV	1.35	CV	1.69
Day-Run	Value	Day-Run	Value
1 - 1	173	1 - 1	236
3 - 1	172	3 - 1	233
4 - 1	174	4 - 1	234
6 - 1	173	6 - 1	232
7 - 1	174	7 - 1	237
8 - 1	176	8 - 1	236
10 - 1	171	10 - 1	235
11 - 1	176	11 - 1	235
13 - 1	178	13 - 1	233
14 - 1	174	14 - 1	230
15 - 1	170	15 - 1	230
16 - 1	173	16 - 1	228
20 - 1	175	20 - 1	240
21 - 1	176	21 - 1	228
22 - 1	177	22 - 1	233
23 - 1	176	23 - 1	233
24 - 1	178	24 - 1	227
28 - 1	174	28 - 1	227
29 - 1	170	29 - 1	230
30 - 1	174	30 - 1	235

Westgard	تفسیر	Triglycerides	آزمایش	اتوآنالایزر BS-800	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	173.19	SD	1.95	Exp. Date	1403/1/17	Lot. No	035

