

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.34	
Mean	34.46	Mean	130.7
SD	1.12	SD	1.75
CV	3.24	CV	1.34
Day-Run	Value	Day-Run	Value
1 - 1	31	1 - 1	130
3 - 1	32.8	3 - 1	129
4 - 1	33.9	4 - 1	133
6 - 1	35.1	6 - 1	131.8
7 - 1	34.4	7 - 1	130.7
8 - 1	35	8 - 1	132.3
10 - 1	36	10 - 1	132.8
11 - 1	34.8	11 - 1	130.4
13 - 1	35	13 - 1	127
14 - 1	35	14 - 1	131
15 - 1	34.9	15 - 1	131.8
16 - 1	35	16 - 1	131
20 - 1	34	20 - 1	130.4
21 - 1	34.4	21 - 1	131.6
22 - 1	33.7	22 - 1	128
23 - 1	36	23 - 1	128
24 - 1	35	24 - 1	133.4
28 - 1	34.8	28 - 1	130.4
29 - 1	34.7	29 - 1	132
30 - 1	33.6	30 - 1	129.4

Westgard	تفسیر	Urea	آزمایش	اتوآنالایزر BS-800	دستگاه	8	ماه	1402	سال
Normal Control	biorex	Mean	34.24	SD	1.16	Exp. Date	1403/1/17	Lot. No	035

