

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	64.1	
3 - 1	70.3	
4 - 1	69.1	
6 - 1	68.9	
7 - 1	64.8	
8 - 1	69.6	
10 - 1	69.6	
11 - 1	64.5	
13 - 1	67.5	
14 - 1	70	
15 - 1	69	
16 - 1	69	
20 - 1	68.1	
21 - 1	67.2	
22 - 1	60	1:3S
23 - 1	62.6	
24 - 1	62.9	
28 - 1	67.6	
29 - 1	68.6	
30 - 1	66	

Westgard	تفسیر	Lipase	آزمایش	BS-800 اتوآنالایزر	دستگاه	8	ماه	1402	سال
Control	biorex	Mean	66.92	SD	2.23	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	67.34	SD	2.45	CV	3.64			

