

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	28.12
SD	1.01
CV	3.61
Day-Run	Value
1 - 1	28.5
2 - 1	27.1
5 - 1	28.6
6 - 1	28.1
8 - 1	27.3
9 - 1	27.5
11 - 1	28.1
12 - 1	30
13 - 1	27.1
14 - 1	27.8
15 - 1	28.2
16 - 1	28
18 - 1	26.8
19 - 1	27
20 - 1	28
22 - 1	27
23 - 1	30
25 - 1	29.7
27 - 1	27.3
28 - 1	29.4
29 - 1	29

High - CV : 2.17	
Mean	137.32
SD	2.98
CV	2.17
Day-Run	Value
1 - 1	132.6
2 - 1	138.6
5 - 1	134.8
6 - 1	136.1
8 - 1	140
9 - 1	141
12 - 1	139.5
13 - 1	133
13 - 2	140
14 - 1	135.1
15 - 1	136.4
16 - 1	133
18 - 1	140
19 - 1	143
20 - 1	140
21 - 1	139
22 - 1	135
23 - 1	140
25 - 1	135.4
27 - 1	133.7
28 - 1	136.9
29 - 1	138

Westgard	تفسیر	SGPT (ALT)	آزمایش	اتوآنالایزر BS-800	دستگاه	9	ماه	1402	سال
Normal Control	biorex	Mean	28.86	SD	1.05	Exp. Date	1403/1/17	Lot. No	035

