

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	37.85
SD	0.73
CV	1.92
Day-Run	Value
1 - 1	38.5
2 - 1	36.8
5 - 1	39
6 - 1	37.5
8 - 1	36.5
11 - 1	39
12 - 1	37.8
13 - 1	37.2
14 - 1	38.6
15 - 1	38.3
16 - 1	37
18 - 1	37.3
19 - 1	38
21 - 1	39
22 - 1	37.6
23 - 1	38
25 - 1	37.8
27 - 1	37.7
28 - 1	37.4
29 - 1	38

High - CV : 1.44	
Mean	145.14
SD	2.09
CV	1.44
Day-Run	Value
1 - 1	142.8
2 - 1	147.5
5 - 1	143
6 - 1	142
11 - 1	146.8
11 - 2	143.1
12 - 1	146.8
13 - 1	143.5
14 - 1	147
15 - 1	143
16 - 1	143.1
18 - 1	146
19 - 1	148
20 - 1	147
21 - 1	148
22 - 1	147
23 - 1	145
25 - 1	143.9
27 - 1	144
28 - 1	147.5
29 - 1	143

Westgard	تفسیر	SGOT (AST)	آزمایش	اتوآنالایزر BS-800	دستگاه	9	ماه	1402	سال
Normal Control	biorex	Mean	37.94	SD	1	Exp. Date	1403/1/17	Lot. No	035

