

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	0.47
SD	0.04
CV	8.43
Day-Run	Value
1 - 1	0.42
2 - 1	0.5
5 - 1	0.44
11 - 1	0.46
12 - 1	0.47
13 - 1	0.51
14 - 1	0.42
15 - 1	0.52
16 - 1	0.44
23 - 1	0.4
25 - 1	0.47
27 - 1	0.48
28 - 1	0.5
29 - 1	0.52

High - CV : 3.44	
Mean	2.25
SD	0.08
CV	3.44
Day-Run	Value
1 - 1	2.31
2 - 1	2.24
11 - 1	2.23
12 - 1	2.3
13 - 1	2.33
14 - 1	2.1
15 - 1	2.33
16 - 1	2.15
23 - 1	2.16
25 - 1	2.27
27 - 1	2.2
28 - 1	2.3
29 - 1	2.33

Westgard	تفسیر	D.Bili	آزمایش	BS-800 اتوآنالایزر	دستگاه	9	ماه	1402	سال
Normal Control	biorex	Mean	0.45	SD	0.05	Exp. Date	1403/1/17	Lot. No	035

