

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

| Normal | |
|---------|-------|
| Mean | 1.01 |
| SD | 0.07 |
| CV | 7.18 |
| Day-Run | Value |
| 1 - 1 | 1 |
| 2 - 1 | 0.95 |
| 11 - 1 | 1 |
| 12 - 1 | 1.1 |
| 13 - 1 | 0.93 |
| 13 - 1 | 0.93 |
| 14 - 1 | 1.05 |
| 15 - 1 | 1.11 |
| 16 - 1 | 0.93 |
| 23 - 1 | 1.1 |
| 25 - 1 | 0.96 |
| 27 - 1 | 1.03 |
| 28 - 1 | 0.93 |
| 29 - 1 | 1.1 |

| High - CV : 2.76 | |
|------------------|-------|
| Mean | 4.27 |
| SD | 0.12 |
| CV | 2.76 |
| Day-Run | Value |
| 1 - 1 | 4.2 |
| 2 - 1 | 4.23 |
| 11 - 1 | 4.18 |
| 12 - 1 | 4.43 |
| 13 - 1 | 4.11 |
| 14 - 1 | 4.49 |
| 15 - 1 | 4.18 |
| 16 - 1 | 4.3 |
| 23 - 1 | 4.3 |
| 25 - 1 | 4.3 |
| 27 - 1 | 4.2 |
| 28 - 1 | 4.4 |
| 29 - 1 | 4.13 |

| | | | | | | | | | |
|----------------|--------|--------|--------|--------------------|--------|-----------|-----------|---------|-----|
| Westgard | تفسیر | Bili T | آزمایش | اتوآنالایزر BS-800 | دستگاه | 9 | ماه | 1402 | سال |
| Normal Control | biorex | Mean | 1.02 | SD | 0.07 | Exp. Date | 1403/1/17 | Lot. No | 035 |

