

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	91.6	
2 - 1	92	
4 - 1	91.3	
5 - 1	90.2	
8 - 1	91.1	
9 - 1	91	
11 - 1	91.7	
12 - 1	92.1	
13 - 1	90.7	
14 - 1	90.3	
15 - 1	91.7	
16 - 1	91	
18 - 1	90.1	
19 - 1	90.3	
20 - 1	91.5	
21 - 1	91.8	
22 - 1	90.2	
23 - 1	91	
25 - 1	91.7	
27 - 1	90.8	
28 - 1	91.03	
29 - 1	90.2	
30 - 1	91.9	

Westgard	تفسیر	C3	آزمایش	BS-480 اتوآنالایزر	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	91.06	SD	0.52	Exp. Date	1402/11/7	Lot. No	641-01-030/2
Next Month	Mean	91.1	SD	0.66	CV	0.72			

