

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	18.2	
2 - 1	17.6	
4 - 1	18	
5 - 1	18.1	
8 - 1	17.5	1:2S
9 - 1	17.8	
11 - 1	18.2	
12 - 1	17.8	
13 - 1	18.31	
14 - 1	18.2	
15 - 1	17.6	
16 - 1	17.8	
18 - 1	17.5	1:2S
19 - 1	18.3	
20 - 1	18.2	
21 - 1	17.7	
22 - 1	17.9	
23 - 1	18.1	
25 - 1	17.5	1:2S
27 - 1	17.9	
28 - 1	18.3	
29 - 1	17.93	
30 - 1	18.3	

Westgard	تفسیر	C4	آزمایش	BS-480 اتوآنالایزر	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	17.94	SD	0.19	Exp. Date	1402/9/14	Lot. No	641-01-030/2
Next Month	Mean	17.95	SD	0.28	CV	1.59			

