

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	49.3	
4 - 1	53.9	
5 - 1	51.1	
6 - 1	54.2	
7 - 1	57.2	1:2S
8 - 1	52.5	
9 - 1	53.1	
11 - 1	57	1:2S
12 - 1	49	
13 - 1	51.6	
14 - 1	49.3	
15 - 1	52	
16 - 1	54.9	
18 - 1	53	
19 - 1	47.8	
20 - 1	54.5	
21 - 1	48.1	
22 - 1	52.8	
27 - 1	48.5	
28 - 1	51.3	
29 - 1	51.5	
30 - 1	52.1	

Westgard	تفسیر	CA 125	آزمایش	2000 ایمولایت XPI-	دستگاه	9	ماه	1402	سال
Control	Seronorm	Mean	50.3	SD	2.34	Exp. Date	1403/10/1	Lot. No	2011909
Next Month	Mean	52.03	SD	2.66	CV	5.11			

