

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	174.95
SD	2.68
CV	1.53
Day-Run	Value
1 - 1	171
2 - 1	171
5 - 1	170
6 - 1	172
7 - 1	173
11 - 1	175
12 - 1	177
13 - 1	178
14 - 1	177
15 - 1	172
16 - 1	176
18 - 1	178
19 - 1	177
20 - 1	175
21 - 1	178
22 - 1	177
25 - 1	176
27 - 1	173
28 - 1	178
29 - 1	175

High - CV : 0.85	
Mean	261.85
SD	2.23
CV	0.85
Day-Run	Value
1 - 1	260
2 - 1	263
5 - 1	262
6 - 1	264
7 - 1	259
11 - 1	258
12 - 1	260
13 - 1	262
14 - 1	262
15 - 1	259
16 - 1	260
18 - 1	263
19 - 1	260
20 - 1	263
21 - 1	266
22 - 1	263
23 - 1	270
25 - 1	265
27 - 1	260
28 - 1	264
29 - 1	264

Westgard	تفسیر	Cholesterol	آزمایش	اتوآنالایزر BS-800	دستگاه	9	ماه	1402	سال
Normal Control	biorex	Mean	173.35	SD	2.58	Exp. Date	1403/1/17	Lot. No	035

