

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.09	
Mean	126.05	Mean	218.95
SD	1.77	SD	4.58
CV	1.41	CV	2.09
Day-Run	Value	Day-Run	Value
1 - 1	126	1 - 1	224
2 - 1	126	2 - 1	221
5 - 1	124	5 - 1	213
6 - 1	124	6 - 1	217
8 - 1	127	11 - 1	230
11 - 1	124	11 - 2	219
12 - 1	128	12 - 1	226
13 - 1	128	13 - 1	218
14 - 1	126	14 - 1	219
15 - 1	126	14 - 2	219
16 - 1	126	15 - 1	224
18 - 1	123	16 - 1	214
19 - 1	123	18 - 1	214
20 - 1	125	20 - 1	215
22 - 1	126	21 - 1	220
23 - 1	127	22 - 1	221
25 - 1	127	23 - 1	218
27 - 1	127	25 - 1	214
28 - 1	126	27 - 1	214
29 - 1	130	28 - 1	215
29 - 2	128	29 - 1	223

Westgard	تفسیر	CPK	آزمایش	اتوآنالایزر BS-800	دستگاه	9	ماه	1402	سال
Normal Control	biorex	Mean	125.55	SD	1.61	Exp. Date	1403/1/17	Lot. No	035

