

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	48.1	
2 - 1	47.3	
4 - 1	48.3	
5 - 1	47.3	
8 - 1	47	
9 - 1	48.1	
11 - 1	47.2	
12 - 1	47.1	
13 - 1	48.3	
14 - 1	46.2	1:2S
15 - 1	47.3	
16 - 1	46.8	
18 - 1	47.5	
19 - 1	48.2	
20 - 1	47.1	
21 - 1	46.9	
22 - 1	48.3	
23 - 1	47.5	
25 - 1	46.8	1:2S
27 - 1	47.6	
28 - 1	48	
29 - 1	46.3	1:2S
30 - 1	47.4	

Westgard	تفسیر	CRP-quantitative	آزمایش	BS-480 اتوآنالایزر	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	47.45	SD	0.54	Exp. Date	1403/6/8	Lot. No	326-2-034/3
Next Month	Mean	47.42	SD	0.62	CV	1.31			

