

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	30.7	
2 - 1	29.8	
4 - 1	30.2	
5 - 1	29.3	
8 - 1	30.1	
9 - 1	31.42	
11 - 1	31.3	
12 - 1	30.7	
13 - 1	30.1	4:1S
14 - 1	29.2	
15 - 1	30.5	
16 - 1	31.2	
18 - 1	30.1	
19 - 1	30.2	1:2S
20 - 1	29.7	
21 - 1	31.1	
22 - 1	30.7	
23 - 1	29.5	
25 - 1	31.1	
27 - 1	31.2	
28 - 1	30.5	
29 - 1	28.6	
30 - 1	30.4	

Westgard	تفسیر	Ceruloplasmin	آزمایش	BS-480 اتوآنالایزر	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	29.85	SD	0.79	Exp. Date	1406/9/14	Lot. No	641-01-030/2
Next Month	Mean	30.37	SD	0.76	CV	2.49			

