

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
11 - 1	24.9	
12 - 1	24.2	
13 - 1	24.8	
14 - 1	25.3	
15 - 1	24.6	
16 - 1	24.7	
18 - 1	25.1	
19 - 1	24.9	
20 - 1	24.78	
21 - 1	25.1	
22 - 1	25.1	
23 - 1	25.2	
25 - 1	24.3	
27 - 1	25	
28 - 1	24.2	
30 - 1	24.6	

Westgard	تفسیر	Cu,urine	آزمایش	Atomic Absorption Spectrometry	دستگاه	9	ماه	1402	سال
Control	آزمایشگاه	Mean	24.75	SD	0.76	Exp. Date	1402/11/30	Lot. No	
Next Month	Mean	24.8	SD	0.35	CV	1.4			

