

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	23.6	
4 - 1	24.7	
5 - 1	24.7	
6 - 1	24.9	
7 - 1	23.8	
8 - 1	24.7	
9 - 1	24.1	
11 - 1	25.5	
12 - 1	24.3	
13 - 1	22.8	1:2S
14 - 1	24.1	
15 - 1	24.4	
16 - 1	25	
18 - 1	24	
20 - 1	23.9	
21 - 1	23.1	
22 - 1	23.4	
25 - 1	25	
27 - 1	23.7	
28 - 1	22.7	1:2S
30 - 1	25.4	

Westgard	تفسیر	FSH	آزمایش	2000 ایمولایت XPI-	دستگاه	9	ماه	1402	سال
Control	آزمایشگاه	Mean	24.75	SD	0.94	Exp. Date	1403/10/1	Lot. No	
Next Month	Mean	24.18	SD	0.8	CV	3.29			

