

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	66.9	
2 - 1	68.9	
5 - 1	70.8	
6 - 1	69.9	
7 - 1	68.6	
8 - 1	65.9	
11 - 1	66.4	
12 - 1	67	
13 - 1	65	1:2S
14 - 1	71.2	
15 - 1	68.9	
16 - 1	67.5	
16 - 1	67.5	
18 - 1	64.7	1:2S
19 - 1	67.2	1:2S
20 - 1	65	1:2S
21 - 1	66	
22 - 1	68	
23 - 1	67.2	
25 - 1	66.7	
27 - 1	70.3	
28 - 1	69	
29 - 1	68.8	

Westgard	تفسیر	Gamma GT	آزمایش	BS-800 اتوآنالایزر	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	68.52	SD	1.64	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	67.71	SD	1.84	CV	2.71			

